



**CoughMonitor
Suite**




Participant Brochure

What is the CoughMonitor Suite?

The CoughMonitor Suite is a tool to monitor cough. It consists of the CoughMonitor watch and the CoughMonitor Companion mobile app. The CoughMonitor watch uses an advanced algorithm to distinguish cough from other sounds without storing audio, which protects your privacy. The Companion app is used to log into the CoughMonitor watch and send data to Hyfe’s Cloud.

CoughMonitor watch components



1.	Round button	Quick press to wake up the screen when it is off
2.	Microphone	Used to capture sounds
3.	Charger port	Used to connect the charging cable
4.	USB charger cable	Used for charging the CoughMonitor watch
5.	Package box	Needed to return the device at the end of the study


CoughMonitor watch screen

1	Battery level
2	Indicates charging in progress
3	Indicates Bluetooth is active
4	Software version
5	Date
6	Time (24hours)
7	Indicates monitoring status
8	CoughMonitor unique identifier



How should I use the CoughMonitor watch?

Every day for 9 consecutive weeks, place the CoughMonitor watch on your wrist and carry out your normal activities as usual.



Whenever possible, avoid covering the CoughMonitor watch with thick clothing.


You should not submerge it in water for more than 30 minutes or deeper than 3 feet

You should not wear the watch below -20°C (4°F) or above 50°C (122°F)

Open the app on your smartphone once a day to check data has synced


Do I need to charge the watch?

Yes, before going to sleep, take off the CoughMonitor watch and leave it charging at your bedside, no further than 3 feet/90cm away from you.



What happens if I forget to charge the CoughMonitor?

If you forget to charge the CoughMonitor watch and the battery drains, take it off and connect it to the charger cable. Make sure Bluetooth is enabled on your phone, and open the Companion app. By simply opening the Companion app, the CoughMonitor watch will restart cough monitoring. The Companion app will show the CoughMonitor as "Connected."



If the Companion app shows "Disconnected," tap the "Try Again" button to force reconnection. Close the Companion app (see instructions in the next page), then open the app again. If after 2 minutes the app still shows the CoughMonitor as "Disconnected," contact your research team.

What should I do if I forget to wear my watch?

As soon as you remember, please put on your watch to wear each day, and charge at night. If you forget to put on your watch, please notify your research team.

If you have any questions, please contact us at:

[INSERT SITE CONTACT INFORMATION]

PHONE NUMBER

EMAIL

ADDRESS


[INSERT CONCIERGE CONTACT INFORMATION]

[INSERT SITE CONTACT INFORMATION]

PHONE NUMBER

EMAIL

ADDRESS



How do I return my device?

Once you have completed the study, you will be asked to open the app on your phone and ensure all data has been synced. You can then put the watch and its accessories back in the original box, place the return label and ship at your nearest convenience.

How do I close the Companion app to troubleshoot if the CoughMonitor shows as 'disconnected'?

Android: Swipe up from the bottom of the screen, hold, then let go. Alternative, press the three-button navigation to switch apps. Swipe right or left to find the app that you want to close. Swipe up on the app.

iPhone: From the Home Screen, swipe up from the bottom of the screen and pause in the middle of the screen. Swipe right or left to find the app that you want to close. Swipe up on the app's preview to close the app.